

Product Datasheet

ACSS1 (Acetyl-K642) antibody (orb1415994)

Description	Rabbit polyclonal antibody to ACSS1 (Acetyl-K642).
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q9NUB1
Dilution Range	WB: Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Expiration Date	12 months from date of receipt.