

Product Datasheet

HSP70 (Acetyl Lys246) antibody (orb1416017)

Description	Rabbit polyclonal antibody to HSP70 (Acetyl Lys246).
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P08107
Dilution Range	WB: Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Expiration Date	12 months from date of receipt.