

Product Datasheet

Tubulin α (Acetyl Lys40) antibody (orb1416039)

Description	Rabbit polyclonal antibody to Tubulin α (Acetyl Lys40).
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q71U36
Dilution Range	WB: Western Blot: 1/500 - 1/2000. IHC-p: 1:100-300 ELISA: 1/20000. Not yet tested in other applications.
Expiration Date	12 months from date of receipt.