

Product Datasheet

Ref-1 (Acetyl Lys6) antibody (orb1416056)

Description	Rabbit polyclonal antibody to Ref-1 (Acetyl Lys6).
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P27695
Dilution Range	WB: Western Blot: 1/500 - 1/2000. IHC-p: 1:100-300 ELISA: 1/20000. Not yet tested in other applications.
Expiration Date	12 months from date of receipt.