

## Product Datasheet

### Histone H3 (Acetyl Lys27) antibody (orb1416065)

|                            |                                                                                                                                                                          |
|----------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to Histone H3 (Acetyl Lys27).                                                                                                                 |
| <b>Species/Host</b>        | Rabbit                                                                                                                                                                   |
| <b>Reactivity</b>          | Human, Mouse, Rat                                                                                                                                                        |
| <b>Conjugation</b>         | Unconjugated                                                                                                                                                             |
| <b>Tested Applications</b> | IF, IHC-P, WB                                                                                                                                                            |
| <b>Note</b>                | For research use only                                                                                                                                                    |
| <b>Clonality</b>           | Polyclonal                                                                                                                                                               |
| <b>Uniprot ID</b>          | <b>Q71DI3</b> , <b>P84243</b> , <b>P68431</b>                                                                                                                            |
| <b>Dilution Range</b>      | WB: Western Blot: 1/500 - 1/2000.<br>Immunohistochemistry: 1/100 - 1/300.<br>Immunofluorescence: 1/200 - 1/1000.<br>ELISA: 1/5000. Not yet tested in other applications. |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                                                                          |