

Product Datasheet

Histone H3 (Acetyl Lys18) antibody (orb1416067)

Description	Rabbit polyclonal antibody to Histone H3 (Acetyl Lys18).
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC-P, WB
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q71DI3 , P84243 , P68431
Dilution Range	WB: Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/5000. Not yet tested in other applications.
Expiration Date	12 months from date of receipt.