

Product Datasheet

Histone H3 (Acetyl Lys9) antibody (orb1416069)

Description	Rabbit polyclonal antibody to Histone H3 (Acetyl Lys9).
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC-P, WB
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q71DI3 , P84243 , P68431
Dilution Range	WB: WB 1:500-2000, IHC-p 1:50-300, IF 1:50-300
Expiration Date	12 months from date of receipt.