

Product Datasheet

Cleaved-PARP-1 (D214) antibody (orb1416092)

Description	Rabbit polyclonal antibody to Cleaved-PARP-1 (D214).
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IF, IHC-P, WB
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P09874
Dilution Range	WB: WB 1:500-2000, IF 1:50-300, IHC 1:50-300
Expiration Date	12 months from date of receipt.