

Product Datasheet

Cleaved-PAR-4 (G48) antibody (orb1416114)

Description	Rabbit polyclonal antibody to Cleaved-PAR-4 (G48).
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q96RI0
Dilution Range	WB: Western Blot: 1/500 - 1/2000. ELISA: 1/5000. Not yet tested in other applications.
Expiration Date	12 months from date of receipt.