

Product Datasheet

Cleaved-Caspase-9 p35 (D315) antibody (orb1416161)

Description	Rabbit polyclonal antibody to Cleaved-Caspase-9 p35 (D315).
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IHC-P, WB
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P55211
Dilution Range	WB: WB 1:500-2000, IHC-p 1:50-300, IF 1:50-300
Expiration Date	12 months from date of receipt.