

Product Datasheet

Cleaved-Caspase-3 p17 (D175) antibody (orb1416168)

Description	Rabbit polyclonal antibody to Cleaved-Caspase-3 p17 (D175).
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC-P, WB
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P42574
Dilution Range	WB: WB 1:500-2000, IHC-p 1:50-300, IF 1:50-300
Expiration Date	12 months from date of receipt.