

Product Datasheet

Cleaved-Caspase-3 p17 (D175) antibody (orb1416169)

Description	Rabbit polyclonal antibody to Cleaved-Caspase-3 p17 (D175).
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P42574
Dilution Range	IHC: Immunohistochemistry: 1/100 - 1/300. ELISA: 1/20000. Not yet tested in other applications.
Expiration Date	12 months from date of receipt.