

Product Datasheet

iNOS antibody (FITC) (orb14179)

Description	iNOS antibody (FITC)
Species/Host	Rabbit
Conjugation	FITC
Tested Applications	ICC, IF
Immunogen	Synthetic peptide (please contact us for details)
Target	iNOS (FITC)
Form/Appearance	Liquid: PBS, 2 mM Sodium Azide, 30% glycerol
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
MW	130 kDa
Uniprot ID	P35228, Q06518, P29477
NCBI	NM_000625.41, NP_000616.31, NP_000616.3, NM_000625.4
Entrez	4843
Dilution Range	IF/ICC:1:50-200 (based on 0.5 mg/ml)
Expiration Date	12 months from date of receipt.