

## Product Datasheet

### ILTV glycoprotein E antibody (orb1419249)

|                            |                                                    |
|----------------------------|----------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to ILTV glycoprotein E. |
| <b>Species/Host</b>        | Rabbit                                             |
| <b>Conjugation</b>         | Unconjugated                                       |
| <b>Tested Applications</b> | IF, IHC-Fr, IHC-P, WB                              |
| <b>Note</b>                | For research use only                              |
| <b>Clonality</b>           | Polyclonal                                         |
| <b>Uniprot ID</b>          | <b>Q0PEE1</b>                                      |
| <b>Expiration Date</b>     | 12 months from date of receipt.                    |