

## Product Datasheet

### C20orf152 antibody (orb1419838)

|                            |                                          |
|----------------------------|------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to C20orf152. |
| <b>Species/Host</b>        | Rabbit                                   |
| <b>Reactivity</b>          | Bovine, Canine, Rat                      |
| <b>Conjugation</b>         | Unconjugated                             |
| <b>Tested Applications</b> | ICC, IF, IHC-Fr, IHC-P, WB               |
| <b>Note</b>                | For research use only                    |
| <b>Clonality</b>           | Polyclonal                               |
| <b>Uniprot ID</b>          | <b>Q96M20</b>                            |
| <b>Expiration Date</b>     | 12 months from date of receipt.          |