

## Product Datasheet

### C1orf69 antibody (orb1419854)

|                            |                                        |
|----------------------------|----------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to C1orf69. |
| <b>Species/Host</b>        | Rabbit                                 |
| <b>Reactivity</b>          | Porcine, Rat                           |
| <b>Conjugation</b>         | Unconjugated                           |
| <b>Tested Applications</b> | ICC, IF, IHC-Fr, IHC-P, WB             |
| <b>Note</b>                | For research use only                  |
| <b>Clonality</b>           | Polyclonal                             |
| <b>Uniprot ID</b>          | <b>Q5T440</b>                          |
| <b>Expiration Date</b>     | 12 months from date of receipt.        |