

Product Datasheet

Alpha-actin-1 antibody (orb1421525)

Description	Mouse monoclonal antibody to Alpha-actin-1.
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IH, IP, WB
Note	For research use only
Clonality	Monoclonal
Uniprot ID	P68133
Dilution Range	WB: WB (1/5000 - 1/10000), IH (1/200 - 1/500), IP (1/100 - 1/200), IH: WB (1/5000 - 1/10000), IH (1/200 - 1/500), IP (1/100 - 1/200), IP: WB (1/5000 - 1/10000), IH (1/200 - 1/500), IP (1/100 - 1/200)
Expiration Date	12 months from date of receipt.