

Product Datasheet

DARPP32 antibody (orb1421635)

Description	Rabbit polyclonal antibody to DARPP32.
Species/Host	Rabbit
Reactivity	Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q9UD71
Dilution Range	WB: WB (1/500 - 1/1000), IH (1/50 - 1/200), IH: WB (1/500 - 1/1000), IH (1/50 - 1/200)
Expiration Date	12 months from date of receipt.