

Product Datasheet

DAP Kinase 3 antibody (orb1421800)

Description	Rabbit polyclonal antibody to DAP Kinase 3.
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IH, WB
Note	For research use only
Clonality	Polyclonal
Uniprot ID	O43293
Dilution Range	WB: WB (1/500 - 1/1000), IH (1/50 - 1/200), IH: WB (1/500 - 1/1000), IH (1/50 - 1/200)
Expiration Date	12 months from date of receipt.