

Product Datasheet

BLNK (pY84) antibody (orb1421826)

Description	Rabbit polyclonal antibody to BLNK (pY84).
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IH, WB
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q8WV28
Dilution Range	WB: WB (1/500 - 1/1000), IH (1/50 - 1/200), IH: WB (1/500 - 1/1000), IH (1/50 - 1/200)
Expiration Date	12 months from date of receipt.