

Product Datasheet

DCLK2 antibody (orb1422239)

Description	Rabbit polyclonal antibody to DCLK2.
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, WB
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q8N568
Dilution Range	WB: WB (1/500 - 1/2000), IF/IC (1/50 - 1/100), IF: WB (1/500 - 1/2000), IF/IC (1/50 - 1/100)
Expiration Date	12 months from date of receipt.