

Product Datasheet

AT1 antibody (orb1422278)

Description	Rabbit polyclonal antibody to AT1.
Species/Host	Rabbit
Reactivity	Human, Mouse, Porcine, Rat, Virus
Conjugation	Unconjugated
Tested Applications	IF, IH, WB
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P30556
Dilution Range	WB: WB (1/500 - 1/2000), IH (1/50 - 1/200), IF/IC (1/50 - 1/100), IH: WB (1/500 - 1/2000), IH (1/50 - 1/200), IF/IC (1/50 - 1/100), IF: WB (1/500 - 1/2000), IH (1/50 - 1/200), IF/IC (1/50 - 1/100)
Expiration Date	12 months from date of receipt.