

Product Datasheet

Alpha-synuclein (pY133) antibody (orb1422576)

Description	Rabbit polyclonal antibody to Alpha-synuclein (pY133).
Species/Host	Rabbit
Reactivity	Human, Mouse, Porcine, Rabbit, Rat, Virus
Conjugation	Unconjugated
Tested Applications	IH, WB
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P37840
Dilution Range	WB: WB (1/500 - 1/1000), IH (1/100 - 1/200), IH: WB (1/500 - 1/1000), IH (1/100 - 1/200)
Expiration Date	12 months from date of receipt.