

Product Datasheet

Histone H3 (AcK27) antibody (orb1422777)

Description	Rabbit polyclonal antibody to Histone H3 (AcK27).
Species/Host	Rabbit
Reactivity	Human, Mouse, Porcine, Rabbit, Rat, Zebrafish
Conjugation	Unconjugated
Tested Applications	IF, IH, WB
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q71DI3 , P84243 , P68431
Dilution Range	WB: WB (1/500 - 1/1000), IH (1/100 - 1/200), IF/IC (1/100 - 1/500), IH: WB (1/500 - 1/1000), IH (1/100 - 1/200), IF/IC (1/100 - 1/500), IF: WB (1/500 - 1/1000), IH (1/100 - 1/200), IF/IC (1/100 - 1/500)
Expiration Date	12 months from date of receipt.