

Product Datasheet

NF-kappaB p65 (pS536) antibody (orb1423401)

Description Rabbit polyclonal antibody to NF-kappaB p65 (pS536).

Species/Host Rabbit

Reactivity Human, Mouse, Porcine, Rat

Conjugation Unconjugated

Tested Applications IH, IP, WB

Note For research use only

Clonality Polyclonal

Uniprot ID [Q04206](#)

Dilution Range WB: WB (1/500 - 1/1000), IH (1/100 - 1/200), IP (1/10 - 1/100), IH: WB (1/500 - 1/1000), IH (1/100 - 1/200), IP (1/10 - 1/100), E: WB (1/500 - 1/1000), IH (1/100 - 1/200), IP (1/10 - 1/100)

Expiration Date 12 months from date of receipt.