

Product Datasheet

NF-kappaB p65 (pS276) antibody (orb1423402)

Description	Rabbit polyclonal antibody to NF-kappaB p65 (pS276).
Species/Host	Rabbit
Reactivity	Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	IH, IP, WB
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q04206
Dilution Range	WB: WB (1/500 - 1/1000), IH (1/100 - 1/200), IP (1/10 - 1/100), IH: WB (1/500 - 1/1000), IH (1/100 - 1/200), IP (1/10 - 1/100), E: WB (1/500 - 1/1000), IH (1/100 - 1/200), IP (1/10 - 1/100)
Expiration Date	12 months from date of receipt.