

Product Datasheet

c-Myc (pT58) antibody (orb1423461)

Description	Rabbit polyclonal antibody to c-Myc (pT58).
Species/Host	Rabbit
Reactivity	Human, Mouse, Porcine, Rat, Virus, Zebrafish
Conjugation	Unconjugated
Tested Applications	IH, IP, WB
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P01106
Dilution Range	WB: WB (1/500 - 1/1000), IH (1/100 - 1/200), IP (1/10 - 1/100), IH: WB (1/500 - 1/1000), IH (1/100 - 1/200), IP (1/10 - 1/100)
Expiration Date	12 months from date of receipt.