

Product Datasheet

IGF1 Receptor antibody (orb1423504)

Description	Rabbit polyclonal antibody to IGF1 Receptor.
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat, Zebrafish
Conjugation	Unconjugated
Tested Applications	IH, WB
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P06213, P08069
Dilution Range	WB: WB (1/500 - 1/1000), IH (1/100 - 1/200), IH: WB (1/500 - 1/1000), IH (1/100 - 1/200)
Expiration Date	12 months from date of receipt.