

Product Datasheet

**Smad2 (phospho-Ser425) (FITC) [Out of stock]
(orb14244)**

Description

Rabbit polyclonal antibody against

Species/Host	Rabbit
Conjugation	FITC
Tested Applications	FC, IF
Immunogen	KLH conjugated synthetic peptide derived from the region around the phosphorylation site of S human SMAD2 (phospho-Ser425)
Target	Smad2 (FITC)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid: 0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Concentration	- 100 µg (in 100 µl): 1 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purified by Protein A
MW	58 kDa
Uniprot ID	Q15796
Dilution Range	Flow-Cyt: 1:50-200, IF: 1:50-200
Expiration Date	12 months from date of receipt.