

Product Datasheet

HIV1 p55 p17 antibody (orb1425)

Description	HIV1 p55 p17 antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-Fr, IHC-P
Immunogen	KLH conjugated synthetic peptide derived from HIV1 p55+p17/HIV1 Pr55Gag (401-460/497 aa)
Target	HIV1 p55+p17
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	17/55 kDa
Uniprot ID	P04591
Expiration Date	12 months from date of receipt.