

Product Datasheet

IGF2R antibody (orb1425092)

Description	Rabbit polyclonal antibody to IGF2R.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthesized Peptide
Target	IGF2R
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	274375
Uniprot ID	P11717
Dilution Range	WB: 1:2000
Expiration Date	12 months from date of receipt.