

Product Datasheet

(DANRE) rab11fip4a antibody (orb1425345)

Description	Rabbit polyclonal antibody to (DANRE) rab11fip4a.
Species/Host	Rabbit
Reactivity	Zebrafish
Conjugation	Unconjugated
Tested Applications	WB
Target	rab11fip4a
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	70946
Uniprot ID	Q3LGD4
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.