

Product Datasheet

STX1B antibody (orb1425377)

Description	Rabbit polyclonal antibody to STX1B.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Target	STX1B
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	33245
Uniprot ID	P61266
Dilution Range	WB: 1:2000
Expiration Date	12 months from date of receipt.