

Product Datasheet

CTSD antibody (orb1426654)

Description	Mouse monoclonal antibody to CTSD.
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Recombinant Protein
Target	CTSD
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	44552
Uniprot ID	P07339
NCBI	NP_001900.1
Dilution Range	WB: 1:2000, IHC-P: 1:25
Expiration Date	12 months from date of receipt.