



Product Datasheet

GC antibody (orb1427645)

Description	Rabbit polyclonal antibody to GC.
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	FC, IHC-P, WB
Immunogen	Synthetic Peptide
Target	GBA
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	59716
Uniprot ID	P04062
NCBI	NP_001005742.1, NP_000148.2, NP_001005741.1
Dilution Range	WB: 1:2000, WB: 1:1000, IHC-P-Leica: 1:500, FC: 1:25
Expiration Date	12 months from date of receipt.



Western blot:All
lanes:Anti-GC
Antibody ...



Western blot analysis
of GC Antibody (Ce...



Immunohistochemical
analysis of human
sk...