



Product Datasheet

LTF antibody (orb1427812)


Description	Rabbit polyclonal antibody to LTF.
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	Recombinant Protein
Target	LTF
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	78182
Uniprot ID	P02788
NCBI	NP_001186078.1, NP_002334.2
Dilution Range	WB: 1:1000, WB: 1:1000, WB: 1:4000, WB: 1:1000-2000, WB: 1:1000, FC: 1:10~50
Expiration Date	12 months from date of receipt.



Western blot analysis of LTF Antibody in...



Western blot analysis of LTF Antibody in...



Western blot:All lanes:Anti-LTF Antibody...