

Product Datasheet

MOSC2 antibody (orb1427824)

Description	Rabbit polyclonal antibody to MOSC2.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC-P, WB
Immunogen	Synthetic Peptide
Target	MTARC2
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	38023
Uniprot ID	Q969Z3
NCBI	NP_060368.2
Dilution Range	WB: 1:1000, WB: 1:1000, IHC-P: 1:50~100, FC: 1:10~50
Expiration Date	12 months from date of receipt.