



Product Datasheet

PFKL (L684) antibody (orb1428100)

Description	Rabbit polyclonal antibody to PFKL (L684).
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Synthetic Peptide
Target	PFKL
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	85018
Uniprot ID	P17858
NCBI	NP_002617.3
Dilution Range	WB: 1:1000, IHC-P: 1:50~100
Expiration Date	12 months from date of receipt.



Western blot analysis of PFKL in HepG2 c...



Immunohistochemical analysis of human ca...