



Product Datasheet

Natriuretic Peptide Receptor C antibody (orb1428125)


Description	Rabbit polyclonal antibody to Natriuretic Peptide Receptor C.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IHC-P, WB
Immunogen	Synthetic Peptide
Target	NPR3
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	59808
Uniprot ID	P17342
NCBI	NP_001191304.1 , NP_001191305.1 , NP_000899.1
Dilution Range	IF: 1:200, WB: 1:1000, IHC-P: 1:10~50
Expiration Date	12 months from date of receipt.



Fluorescent analysis of HeLa cells stain...



Western blot:All lanes:Anti-ANPC Antibod...



Immunohistochemical analysis of human ki...