



Product Datasheet

ROR1 antibody (orb1428592)


Description	Rabbit polyclonal antibody to ROR1.
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	FC, IHC-P, WB
Immunogen	Recombinant Protein
Target	ROR1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	104283
Uniprot ID	Q01973
NCBI	NP_001077061.1 , NP_005003.2
Dilution Range	WB: 1:4000, WB: 1:4000, IHC-P: 1:100, IHC-P: 1:100, FC: 1:25
Expiration Date	12 months from date of receipt.



Western blot:All
lanes:Anti-ROR1
Antibod...



Western blot:Anti-
ROR1 Antibody at
1:400...



Immunohistochemical
analysis of
orb14285...