



Product Datasheet

EphA2 antibody (orb1428723)


Description	Rabbit polyclonal antibody to EphA2.
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	FC, IHC-P, WB
Immunogen	Synthetic Peptide
Target	EPHA2
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	108266
Uniprot ID	P29317
NCBI	NP_004422.2
Dilution Range	WB: 1:1000, WB: 1:1000, WB: 1:1000, WB: 1:1000, IHC-P: 1:50~100, FC: 1:10~50
Expiration Date	12 months from date of receipt.



Western blot: Anti-EPHA2 Antibody (T45) a...



Western blot analysis of hEPHA2-T45 in M...



Western blot analysis of hEPHA2-T45 in m...