

Product Datasheet

CATB antibody (orb1428974)

Description	Rabbit polyclonal antibody to CATB.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC-P, WB
Immunogen	Synthetic Peptide
Target	CTSB
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	37822
Uniprot ID	P07858
NCBI	NP_680090.1, NP_001899.1, NP_680093.1, NP_680092.1, NP_680091.1
Dilution Range	WB: 1:1000, IHC-P: 1:10~50, FC: 1:10~50
Expiration Date	12 months from date of receipt.