



Product Datasheet

PDK4 antibody (orb1429323)


Description	Rabbit polyclonal antibody to PDK4.
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IF, IHC-P, WB
Immunogen	Synthetic Peptide
Target	PDK4
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	46469
Uniprot ID	Q16654
NCBI	NP_002603.1
Dilution Range	IF: 1:100, WB: 1:2000, WB: 1:2000, WB: 1:2000, IHC-P: 1:10~50, IHC-P: 1:25
Expiration Date	12 months from date of receipt.



Fluorescent
confocal
image of
HeLa cells...



Western
blot:Anti-
PDK4
Antibody (C-
term)...



Western
blot:All
lanes:Anti-
PDK4
Antibod...