



Product Datasheet

CD14 antibody (orb1429922)


Description	Rabbit polyclonal antibody to CD14.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC-P, WB
Immunogen	Synthetic Peptide
Target	CD14
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	40076
Uniprot ID	P08571
NCBI	NP_000582.1 , NP_001167575.1 , NP_001167576.1 , NP_001035110.1
Dilution Range	IF: 1:25, WB: 1:2000, WB: 1:1000, WB: 1:1000, WB: 1:2000, IHC-P-Leica: 1:1000, FC: 1:25
Expiration Date	12 months from date of receipt.



Immunofluorescent analysis of A549 cells...



Western blot:All lanes:Anti-CD14 Antibod...



Western blot:All lanes:Anti-CD14 Antibod...