

Product Datasheet

SCDGFB antibody (orb1429967)

Description	Rabbit polyclonal antibody to SCDGFB.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Synthetic Peptide
Target	Pdgfd
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	42809
Uniprot ID	Q9EQT1
NCBI	NP_076452.1
Dilution Range	WB: 1:1000, WB: 1:1000, IHC-P: 1:10~50, IHC-P: 1:10~50
Expiration Date	12 months from date of receipt.