



Product Datasheet

WISP2 antibody (orb1429974)

Description	Rabbit polyclonal antibody to WISP2.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Synthetic Peptide
Target	CCN5
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	26825
Uniprot ID	O76076
NCBI	NP_003872.1
Dilution Range	WB: 1:1000, IHC-P: 1:50~100
Expiration Date	12 months from date of receipt.



Western blot analysis of anti-WISP2 Anti...



Immunohistochemical analysis of human sk...