

Product Datasheet

Presenilin 1 (PSEN1) antibody (orb1429991)

Description	Rabbit polyclonal antibody to Presenilin 1 (PSEN1).
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC-P, WB
Immunogen	Synthetic Peptide
Target	PSEN1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	52668
Uniprot ID	P49768
NCBI	NP_000012.1, NP_015557.2
Dilution Range	IF: 1:10~50, WB: 1:1000, IHC-P: 1:50~100, FC: 1:10~50
Expiration Date	12 months from date of receipt.