



Product Datasheet

SLC10A1 antibody (orb1430246)

Description	Rabbit polyclonal antibody to SLC10A1.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	Synthetic Peptide
Target	SLC10A1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	38119
Uniprot ID	Q14973
NCBI	NP_003040.1
Dilution Range	WB: 1:1000, FC: 1:25
Expiration Date	12 months from date of receipt.



Western blot: All lanes: Anti-SLC10A1
Anti...



Flow cytometric analysis of HepG2 cells
...