

Product Datasheet

ACOX1 antibody (orb1431550)

Description	Rabbit polyclonal antibody to ACOX1.
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	Synthetic Peptide
Target	ACOX1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	74424
Uniprot ID	Q15067
NCBI	NP_004026.2 , NP_001171968.1 , NP_009223.2
Dilution Range	IHC-P: 1:10~50
Expiration Date	12 months from date of receipt.