



Product Datasheet

Nestin (S1409) antibody (orb1431842)


Description	Rabbit polyclonal antibody to Nestin (S1409).
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC-P, WB
Immunogen	Synthetic Peptide
Target	NES
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	177439
Uniprot ID	P48681
NCBI	NP_006608.1
Dilution Range	IF: 1:200, IF: 1:25, WB: 1:2000, WB: 1:2000, IHC-P: 1:25, IHC-P: 1:25, IHC-P: 1:25
Expiration Date	12 months from date of receipt.



Fluorescent confocal image of SH-SY5Y ce...



Immunofluorescent analysis of U-2 OS (hu...



Western blot: All lanes: Anti-Nestin Antib...